



## Bay Navigator Icon Installer Release Notes



These Release Notes contain important information for all Medtech Users. Please ensure that they are circulated amongst all your staff. We suggest that these notes are filed safely for future reference.

## Table of Contents

<b>Adding the Bay Navigator Icon to a Medtech32 Database .....</b>	<b>3</b>
<b>Configuring the Bay Navigator Icon on users Toolbar .....</b>	<b>3</b>

## Adding the Bay Navigator Icon to a Medtech32 Database

1. Ensure that the Bay Navigator Icon Update.exe is available
2. Place the Bay Navigator Icon Update.exe onto the Medtech Server Machine
3. Double Click on the Bay Navigator Icon Update.exe to initiate the installation.
4. The user will be prompted once the installation has completed and clicking on OK will close the Update Complete prompt.


## Configuring the Bay Navigator Icon on users Toolbar

The icons will then be available in the Medtech Toolbar setup to enable the user to select them for display on the Toolbar.



Select the Bay Navigator Icon and move it from the Un Selected to Selected side of the Toolbar setup for the user.

Once completed the icon will be available on the Toolbar when any patient is selected.



When hovering over the icons the Toolbar hints will advise you what URL the icon is linked to.

When clicking on the Bay Navigator Icon on the Toolbar, the selected URL will be launched.

Should you have any questions with regard to the information contained in these Release Notes please don't hesitate to contact Medtech Support on **0800 2 633 832 (Option 1)**.